

## Anti-KCNQ2 Antibody Picoband® PE Conjugated

Catalog Number: A01259-1-PE

### About KCNQ2

KCNQ2, also called kv7.2, is a potassium channel protein coded for by the gene KCNQ2. It is mapped to 20q13.33. The KCNQ2 gene encodes a voltage-gated potassium channel that is expressed in the brain. Expression of human KCNQ2 in *Xenopus laevis* oocytes led to potassium-selective currents that activated slowly with depolarization. Defects in this gene are a cause of benign familial neonatal convulsions type 1 (BFNC), also known as epilepsy, benign neonatal type 1 (EBN1). At least five transcript variants encoding five different isoforms have been found for this gene.

### Overview

Product Name	Anti-KCNQ2 Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O43526

### Technical Details

Immunogen	E.coli-derived human KCNQ2 recombinant protein (Position: M371-D408).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KCNQ2 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.