

Anti-Phospho-KCNQ2/3/4/5 (T217/246/223/251) Antibody

Catalog Number: A01259T217

Overview

Product Name	Anti-Phospho-KCNQ2/3/4/5 (T217/246/223/251) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-KCNQ2/3/4/5 (T217/246/223/251) Antibody catalog # A01259T217. Tested in ELISA, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IHC
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43526

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human Kv7.3/KCNQ3 around the phosphorylation site of Thr246. AA range:191-240
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:100-1:300 ELISA 1:20000

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Phospho-KCNQ2/3/4/5 (T217/246/223/251) Antibody