

## Anti-Mitochondrial Pyruvate dehydrogenase kinase 1/PDK1 Antibody Picoband® FITC Conjugated

Catalog Number: A01268-1-FITC

### About PDK1

PDK1 (Pyruvate Dehydrogenase Kinase Isoenzyme 1), is an enzyme that in humans is encoded by the PDK1 gene. To find human PDKs, Gudi et al. (1995) used oligonucleotide primers to screen a human liver cDNA library by PCR. They identified and reported the deduced amino acid sequences of PDK1, PDK2, and PDK3. The human PDK1 gene encodes a protein with a predicted molecular mass of 49,244 Da that shares 93% identity with that of the rat PDK1 gene. Northern blot analysis showed that the PDK1 message was expressed predominantly in the heart.

### Overview

Product Name	Anti-Mitochondrial Pyruvate dehydrogenase kinase 1/PDK1 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q15118

### Technical Details

Immunogen	E. coli-derived human PDK1 recombinant protein (Position: V107-Y336).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

---

1. PubMed ID: 29845188, lncRNA Ftx promotes aerobic glycolysis and tumor progression through the PPAR $\gamma$  pathway in hepatocellular carcinoma

Visit [bosterbio.com/anti-pdk1-picoband-trade-antibody-a01268-1-boster.html](https://bosterbio.com/anti-pdk1-picoband-trade-antibody-a01268-1-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mitochondrial Pyruvate dehydrogenase kinase 1/PDK1 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.