

Anti-Mitochondrial Pyruvate dehydrogenase kinase 1/PDK1 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A01268-1-Fluoro594

About PDK1

PDK1 (Pyruvate Dehydrogenase Kinase Isoenzyme 1), is an enzyme that in humans is encoded by the PDK1 gene. To find human PDKs, Gudi et al. (1995) used oligonucleotide primers to screen a human liver cDNA library by PCR. They identified and reported the deduced amino acid sequences of PDK1, PDK2, and PDK3. The human PDK1 gene encodes a protein with a predicted molecular mass of 49,244 Da that shares 93% identity with that of the rat PDK1 gene. Northern blot analysis showed that the PDK1 message was expressed predominantly in the heart.

Overview

| | |
|----------------------|--|
| Product Name | Anti-Mitochondrial Pyruvate dehydrogenase kinase 1/PDK1 Antibody Picoband® Fluoro594 Conjugated |
| Reactive Species | Human, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q15118 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E. coli-derived human PDK1 recombinant protein (Position: V107-Y336). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

1 Publications Citing This Product

1. PubMed ID: 29845188, lncRNA Ftx promotes aerobic glycolysis and tumor progression through the PPAR γ pathway in hepatocellular carcinoma

Visit bosterbio.com/anti-pdk1-picoband-trade-antibody-a01268-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mitochondrial Pyruvate dehydrogenase kinase 1/PDK1 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.