

Anti-HNF1 beta (HNF1B) Rabbit Polyclonal Antibody

Catalog Number: A01270


Overview

Product Name	Anti-HNF1 beta (HNF1B) Rabbit Polyclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Rabbit polyclonal HNF1B Antibody. Catalog# A01270. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Rabbit
Uniprot ID	P35680


Technical Details

Immunogen	Synthetic peptide conjugated to KLH derived from within residues 50-150 of Human HNF1B.
Isotype	IgG
Concentration	1 mg/ml
Purification	Purified from the immunized serum by affinity chromatography (Protein A/G)
Suggested Dilutions	WB 1:500


Anti-HNF1 beta (HNF1B) Rabbit Polyclonal Antibody (A01270) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HNF1B (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with Rabbit poly.



Immunohistochemical staining of paraffin-embedded Human pancreas tissue within the normal limits using anti-HNF1B mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)



Immunohistochemical staining of paraffin-embedded Human Kidney tissue within the normal limits using anti-HNF1B mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HNF1 beta (HNF1B) Rabbit Polyclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.