

Anti-Flt-4 Antibody

Catalog Number: A01276-1


Overview

Product Name	Anti-Flt-4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Flt-4 Antibody catalog # A01276-1. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35916


Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 640-700
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-2000 IHC 1:500-200 ELISA 1:10000-20000 IF 1:50-200

Anti-Flt-4 Antibody (A01276-1) Images



Western blot analysis of 293T and mouse kidney lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000.



Immunohistochemical analysis of paraffin-embedded human placenta, antibody was diluted at 1:200.

1 Publications Citing This Product

1. PubMed ID: 24040410, Li C, Fan J, Song X, Zhang B, Chen Y, Li C, Mi K, Ma H, Song Y, Tao X, Li G. Plos One. 2013 Sep 11;8(9):E75388. Doi: 10.1371/Journal.Pone.0075388. Ecollection 2013. Expression Of Angiopoietin-2 And Vascular Endothelial Growth Factor Receptor-3 Cor...

Visit bosterbio.com/anti-flt-4-antibody-a01276-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Flt-4 Antibody

For Research Use Only. Not for use in diagnostic procedures.