

## Anti-ADP-ribosylation factor 1 ARF1 Antibody

Catalog Number: A01279


### Overview

Product Name	Anti-ADP-ribosylation factor 1 ARF1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ADP-ribosylation factor 1 ARF1 Antibody catalog # A01279. Tested in WB,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P84077

### Technical Details

Immunogen	A synthetic peptide of human ARF1.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IF: 1:50-1:200

## Anti-ADP-ribosylation factor 1 ARF1 Antibody (A01279) Images



Western blot analysis of ARF1 expression in HT29 (A), SKOV3 (B), mouse brain (C), rat liver (D) whole cell lysates.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ADP-ribosylation factor 1 ARF1 Antibody

For Research Use Only. Not for use in diagnostic procedures.