

## Anti-SMC5 Antibody Picoband® FITC Conjugated

Catalog Number: A01285-FITC

### About SMC5

Structural maintenance of chromosomes protein 5 is a protein encoded by the SMC5 gene in human. Predicted to enable ATP binding activity. Involved in several processes, including DNA recombination; cellular senescence; and positive regulation of maintenance of mitotic sister chromatid cohesion. Located in cell junction; chromosome; and nuclear body. Part of Smc5-Smc6 complex.

### Overview

Product Name	Anti-SMC5 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8IY18

### Technical Details

Immunogen	E.coli-derived human SMC5 recombinant protein (Position: Q14-Q1089).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SMC5 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.