

## Anti-Cyclophilin A PPIA Antibody

Catalog Number: A01308-1


### Overview

Product Name	Anti-Cyclophilin A PPIA Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cyclophilin A PPIA Antibody catalog # A01308-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P62937


### Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 IHC: 1:50-1:200

## Anti-Cyclophilin A PPIA Antibody (A01308-1) Images



Western blot analysis of Cyclophilin A on different cells lysates using anti-Cyclophilin A antibody at 1/1,000 dilution.  
Positive control: Lane 1: Hela Lane 2: A431 Lane 3: HepG2



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using anti-Cyclophilin A antibody.  
Counter stained with hematoxylin.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cyclophilin A PPIA Antibody

For Research Use Only. Not for use in diagnostic procedures.