

Anti-Protein NDRG1 NDRG1 Antibody

Catalog Number: A01327-2


Overview

Product Name	Anti-Protein NDRG1 NDRG1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Protein NDRG1 NDRG1 Antibody catalog # A01327-2. Tested in WB,ICC/IF,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92597


Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200 IP: 1:50

Anti-Protein NDRG1 NDRG1 Antibody (A01327-2) Images



Western blot analysis of NDRG1 on SW480 cell (1) and mouse eye tissue (2) lysate using anti-NDRG1 antibody at 1/1,000 dilution.



ICC staining NDRG1 in MCF-7 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protein NDRG1 NDRG1 Antibody

For Research Use Only. Not for use in diagnostic procedures.