

Anti-mGluR5 GRM5 Antibody

Catalog Number: A01338

Overview

Product Name	Anti-mGluR5 GRM5 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-mGluR5 GRM5 Antibody catalog # A01338. Tested in IF, IHC, ICC applications. This antibody reacts with Human, Mouse.
Application	IF, IHC, ICC
Clonality	Polyclonal
Storage Instructions	Store at 4°C in the dark
Host	Chicken
Uniprot ID	P41594

Technical Details

Immunogen	Synthetic peptide (keyhole limpet hemocyanin -KLH conjugate) corresponding to different regions of the metabotropic glutamate receptor type 5 gene product, but are shared between the human (NP_001137303.1, NCBI) and mouse (Q3UVX5.2, NCBI) sequences.
Isotype	Chicken IgY
Form	Liquid
Purification	Affinity Purified
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Anti-mGluR5 GRM5 Antibody (A01338) Images




Figure 1. Immunofluorescent staining data of GRM5 using Anti-mGluR5 GRM5 Antibody (A01338).

Adult mouse cochlear ganglion stained with mGluR5 (green) and Tuj-1 (red)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-mGluR5 GRM5 Antibody