

Anti-Furin Antibody

Catalog Number: A01344-1


Overview

Product Name	Anti-Furin Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Furin Antibody catalog # A01344-1. Tested in WB,IHC,ICC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P09958


Technical Details

Immunogen	Recombinant protein within human Furin aa 200-400
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50 ICC: 1:50-1:100

Anti-Furin Antibody (A01344-1) Images



Western blot analysis of Furin on different cell lysates using anti-Furin antibody at 1/500 dilution. Positive control: Lane 1: HepG2 Lane 2: HeLa Lane 3: HeLa Lane 4: MCF-7



Immunohistochemical analysis of paraffin-embedded human liver tissue using anti-Furin antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Furin Antibody

For Research Use Only. Not for use in diagnostic procedures.