

## Anti-Cytokeratin 18 (S48) KRT18 Antibody

Catalog Number: A01357S48


### Overview

Product Name	Anti-Cytokeratin 18 (S48) KRT18 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Cytokeratin 18 (S48) KRT18 Antibody catalog # A01357S48. Tested in WB,IHC,IF applications. This antibody reacts with Human.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P05783


### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 21-70 of Human Cytokeratin 18.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

## Anti-Cytokeratin 18 (S48) KRT18 Antibody (A01357S48) Images



Western blot (WB) analysis of Cytokeratin 18 (S48) polyclonal antibody at 1:500 dilution Lane1:HepG2 whole cell lysate(40ug) Lane2:L02 whole cell lysate(40ug) Lane3:A549 whole cell lysate(40ug) Lane4:PC12 whole cell lysate(40ug) Lane5:MEF whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of Cytokeratin 18 (S48) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytokeratin 18 (S48) KRT18 Antibody

For Research Use Only. Not for use in diagnostic procedures.