

Anti-STAT2 Antibody Picoband® HRP Conjugated

Catalog Number: A01360-HRP

About STAT2

Signal transducer and activator of transcription 2 (STAT2) is a protein that in humans is encoded by the STAT2 gene. The protein encoded by this gene is a member of the STAT protein family. The International Radiation Hybrid Mapping Consortium mapped the STAT2 gene to chromosome 12. STAT2 is a transcription factor critical to the signal transduction pathway of type I interferons. ISGF3 (STAT2) assembly involves p48 functioning as an adaptor protein to recruit Stat1 and Stat2 to an IFN- α -stimulated response element, Stat2 contributes a potent transactivation domain but is unable to directly contact DNA, while Stat1 stabilizes the heteromeric complex by contacting DNA directly.

Overview

Product Name	Anti-STAT2 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P52630

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human STAT2, which shares 73.5% and 76.5% amino acid (aa) sequence identity with mouse and rat STAT2, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-STAT2 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.