

Anti-Dystroglycan DAG1 Antibody

Catalog Number: A01369


Overview

Product Name	Anti-Dystroglycan DAG1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Dystroglycan DAG1 Antibody catalog # A01369. Tested in WB,ICC/IF,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14118


Technical Details

Immunogen	Recombinant fusion protein of human DAG1(NP_004384.4).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 ICC/IF: 1:20-1:100 IP: 1:50-1:200

Anti-Dystroglycan DAG1 Antibody (A01369) Images



Western blot analysis of extracts of various cell lines, using DAG1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunoprecipitation analysis of 300ug extracts of HeLa cells using 3ug DAG1 antibody. Western blot was performed from the immunoprecipitate using DAG1 antibody at a dilution of 1:1000.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Dystroglycan DAG1 Antibody

For Research Use Only. Not for use in diagnostic procedures.