

Anti-VAPB Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A01372-Fluoro550

About VAPB

Vesicle-associated membrane protein-associated protein B/C is a protein that in humans is encoded by the VAPB gene. The VAPB gene is found on the 20th human chromosome. The protein encoded by this gene is a type IV membrane protein found in plasma and intracellular vesicle membranes. The encoded protein is found as a homodimer and as a heterodimer with VAPA. This protein also can interact with VAMP1 and VAMP2 and may be involved in vesicle trafficking.

Overview

| | |
|----------------------|---|
| Product Name | Anti-VAPB Antibody Picoband® Fluoro550 Conjugated |
| Reactive Species | Human, Mouse |
| Application | Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | O95292 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E. coli-derived human VAPB recombinant protein (Position: A2-R55). Human VAPB shares 98.1% and 100% amino acid (aa) sequence identity with mouse and rat VAPB, respectively. |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-VAPB Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.