

Anti-DNA pol iota Antibody

Catalog Number: A01376


Overview

Product Name	Anti-DNA pol iota Antibody
Reactive Species	Human
Description	Boster Bio Anti-DNA pol iota Antibody catalog # A01376. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UNA4

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human POLI. AA range:641-690
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL.
Purification	Immunogen affinity purified
Suggested Dilutions	WB, 1:500-1:2000

Anti-DNA pol iota Antibody (A01376) Images



Western Blot (WB) analysis of HeLa cells using DNA pol iota Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DNA pol iota Antibody