

Anti-Periostin/Postn Antibody Picoband® HRP Conjugated

Catalog Number: A01378-HRP

About Postn

POSTN (Periostin), also known as PN or OSF2, is a protein that in humans is encoded by the POSTN gene. The International Radiation Hybrid Mapping Consortium mapped the POSTN gene to chromosome 13. It was found that purified recombinant PN supported adhesion of ovarian epithelial cells. And periostin was overexpressed by the majority of human primary breast cancers examined. Transfected tumor cell lines overexpressing periostin showed accelerated growth and angiogenesis as xenografts in immunocompromised animals. It was shown that extracellular periostin induced reentry of differentiated mammalian cardiomyocytes into the cell cycle.

Overview

Product Name	Anti-Periostin/Postn Antibody Picoband® HRP Conjugated
Reactive Species	Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	D3ZAF5

Technical Details

Immunogen	E. coli-derived rat Periostin recombinant protein (Position: K577-Q810).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Periostin/Postn Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.