

## Anti-Trefoil factor 1 TFF1 Antibody

Catalog Number: A01391-3

### Overview

Product Name	Anti-Trefoil factor 1 TFF1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Trefoil factor 1 TFF1 Antibody catalog # A01391-3. Tested in WB,ICC/IF,IHC,IP,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P04155

### Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 ICC/IF: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:100 IP: 1:10-1:50

## Anti-Trefoil factor 1 TFF1 Antibody (A01391-3) Images




Figure 1. Immunocytochemistry staining of TFF1 using Anti-Trefoil factor 1 TFF1 Antibody (A01391-3).

ICC staining Estrogen Inducible Protein pS2 in A549 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde




Figure 2. Immunocytochemistry staining of TFF1 using Anti-Trefoil factor 1 TFF1 Antibody (A01391-3).

ICC staining Estrogen Inducible Protein pS2 in MCF-7 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Trefoil factor 1 TFF1 Antibody