

Anti-RAD1 Antibody

Catalog Number: A01396-1


Overview

Product Name	Anti-RAD1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-RAD1 Antibody catalog # A01396-1. Tested in IHC, IF, WB applications. This antibody reacts with Human, Mouse.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O60671

Technical Details

Immunogen	Synthesized peptide derived from human RAD1
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	IHC 1:50-200 WB 1:500-2000 IF 1:50-200

Anti-RAD1 Antibody (A01396-1) Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4°C over night

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAD1 Antibody

For Research Use Only. Not for use in diagnostic procedures.