

## Anti-WNK1 (phospho-T60) Antibody

Catalog Number: A01422T60


### Overview

Product Name	Anti-WNK1 (phospho-T60) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-WNK1 (phospho-T60) Antibody catalog # A01422T60. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H4A3


### Technical Details

Immunogen	Synthetic phosphopeptide derived from human WNK1 around the phosphorylation site of Threonine 60.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-WNK1 (phospho-T60) Antibody (A01422T60) Images



Western blot (WB) analysis of p-WNK1 (T60) pAb at 1:500 dilution Lane1:HEK293T whole cell lysate(40ug)  
Lane2:HEK293T treated with EGF(100ng/ml,30 minutes) whole cell lysate(40ug) Lane3:PC12 whole cell lysate(40ug)  
Lane4:AML-12 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of p-WNK1 (T60) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-WNK1 (phospho-T60) Antibody

For Research Use Only. Not for use in diagnostic procedures.