

Anti-Jagged 2/JAG2 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A01428-1-Fluoro488

About JAG2

Protein jagged-2, also known as SER2, is a protein that in humans is encoded by the JAG2 gene. It is mapped to 14q32.33. The protein encoded by this gene is one of several ligands that activate Notch and related receptors. Notch ligands Delta1 and JAG2 are subject to presenilin-dependent, intramembranous gamma-secretase processing, resulting in the production of soluble intracellular derivatives. JAG2 can promote the growth and dissemination of uveal melanoma cells.

Overview

Product Name	Anti-Jagged 2/JAG2 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y219

Technical Details

Immunogen	E.coli-derived human Jagged 2/JAG2 recombinant protein (Position: Q22-N1153). Human JAG2 shares 91.1% and 91% amino acid (aa) sequence identity with mouse and rat JAG2, respectively.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Jagged 2/JAG2 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.