

## Anti-Cytokeratin 14/KRT14 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A01432-Fluoro647

### About KRT14

Cytokeratin 14, also known as cytokeratin-14 (CK-14) or keratin-14 (KRT14), is a member of the type I keratin family of intermediate filament proteins. In humans it is encoded by the KRT14 gene. This gene product, a type I keratin, is usually found as a heterotetramer with two keratin 5 molecules, a type II keratin. Mutations in the genes for these keratins are associated with epidermolysis bullosa simplex. At least one pseudogene has been identified at 17p12-p11.

### Overview

Product Name	Anti-Cytokeratin 14/KRT14 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IHC, IHC-F, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P02533

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Cytokeratin 14, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

### 3 Publications Citing This Product

---

1. PubMed ID: 10.1042/CBI20090453, Recombinant PBD $\beta$ 1 (porcine beta $\beta$ defensin 1) expressed in the milk by transplanting transgenic mES $\beta$ like $\beta$ derived cells into mouse mammary gland
2. PubMed ID: 10.1007/s11626-008-9092-2, Long-term culture of keratinocyte-like cells derived from mouse embryonic stem cells
3. PubMed ID: 10.1016/j.exer.2006.11.014, Establishment of a corneal epithelial cell line spontaneously derived from human limbal cells

Visit [bosterbio.com/anti-cytokeratin-14-picoband-trade-antibody-a01432-boster.html](http://bosterbio.com/anti-cytokeratin-14-picoband-trade-antibody-a01432-boster.html) to see all 3 publications.

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytokeratin 14/KRT14 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.