

Anti-Annexin A1/ANXA1 Antibody Picoband® Biotin Conjugated

Catalog Number: A01451-1-Biotin

About Annexin A1/ANXA1

Annexin A1, also known as lipocortin I, is a protein that is encoded by the ANXA1 gene in humans. This gene encodes a membrane-localized protein that binds phospholipids. This protein inhibits phospholipase A2 and has anti-inflammatory activity. Loss of function or expression of this gene has been detected in multiple tumors.

Overview

Product Name	Anti-Annexin A1/ANXA1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P04083

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human ANXA1, which shares 82.4% and 76.5% amino acid (aa) sequence identity with mouse and rat ANXA1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Annexin A1/ANXA1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.