

Anti-Cathepsin B (N217) CTSB Antibody

Catalog Number: A01456-2


Overview

| | |
|----------------------|--|
| Product Name | Anti-Cathepsin B (N217) CTSB Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Cathepsin B (N217) CTSB Antibody catalog # A01456-2. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat. |
| Application | IHC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P07858 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Synthetic peptide, corresponding to amino acids 181-230 of Human Cathepsin B. |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000 IHC: 1:50-1:200 |

Anti-Cathepsin B (N217) CTSB Antibody (A01456-2) Images



Western blot (WB) analysis of Cathepsin B (N217) pAb at 1:2000 dilution Lane1:The Embryo tissue lysate of Mouse(40ug) Lane2:The Lung tissue lysate of Rat(40ug) Lane3:HCT116 whole cell lysate(40ug) Lane4:A549 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of Cathepsin B (N217) pAb in paraffin-embedded human brain tissue.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cathepsin B (N217) CTSB Antibody

For Research Use Only. Not for use in diagnostic procedures.