

Anti-GALT Antibody Picoband® Biotin Conjugated

Catalog Number: A01460-1-Biotin

About GALT

Galactose-1-phosphate uridyl transferase (GALT) catalyzes the second step of the Leloir pathway of galactose metabolism, namely the conversion of UDP-glucose + galactose-1-phosphate to glucose-1-phosphate + UDP-galactose. The absence of this enzyme results in classic galactosemia in humans and can be fatal in the newborn period if lactose is not removed from the diet. The pathophysiology of galactosemia has not been clearly defined. Two transcript variants encoding different isoforms have been found for this gene.

Overview

| | |
|----------------------|--|
| Product Name | Anti-GALT Antibody Picoband® Biotin Conjugated |
| Reactive Species | Human, Rat |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | P07902 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E. coli-derived human GALT recombinant protein (Position: Q188-A379). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Biotin |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GALT Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.