

## Anti-Glutathione S-transferase A1 GSTA1 Antibody

Catalog Number: A01462

### Overview


Product Name	Anti-Glutathione S-transferase A1 GSTA1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Glutathione S-transferase A1 GSTA1 Antibody catalog # A01462. Tested in WB, IP, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, IHC, WB
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P08263

### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human GSTA1. AA range:91-140
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ELISA 1:10000-20000 IF 1:50-200

## Anti-Glutathione S-transferase A1 GSTA1 Antibody (A01462) Images

---



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Glutathione S-transferase A1 GSTA1 Antibody

For Research Use Only. Not for use in diagnostic procedures.