

## Anti-GDF-9 Antibody

Catalog Number: A01466


### Overview

Product Name	Anti-GDF-9 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GDF-9 Antibody catalog # A01466. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O60383

### Technical Details

Immunogen	Recombinant fusion protein of human GDF9(NP_005251.1).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

## Anti-GDF-9 Antibody (A01466) Images



Western blot analysis of extracts of SH-SY5Y cells, using GDF9 antibody at 1:2073 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 3s.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-GDF-9 Antibody

For Research Use Only. Not for use in diagnostic procedures.