

## Anti-S6K1/RPS6KB1 Antibody Picoband® APC Conjugated

Catalog Number: A01475-2-APC

### About RPS6KB1

Ribosomal protein S6 kinase beta-1 (S6K1), also known as p70S6 kinase (p70S6K, p70-S6K), is an enzyme (specifically, a protein kinase) that in humans is encoded by the RPS6KB1 gene. This gene encodes a member of the ribosomal S6 kinase family of serine/threonine kinases. The encoded protein responds to mTOR (mammalian target of rapamycin) signaling to promote protein synthesis, cell growth, and cell proliferation. Activity of this gene has been associated with human cancer. Alternatively spliced transcript variants have been observed. The use of alternative translation start sites results in isoforms with longer or shorter N-termini which may differ in their subcellular localizations. There are two pseudogenes for this gene on chromosome 17.

### Overview

Product Name	Anti-S6K1/RPS6KB1 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4, 0.02% NaN3.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P23443

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human S6K1/RPS6KB1, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 10.1007/s10495-017-1391-1, Multi-targeted therapy of everolimus in Kaposi's sarcoma associated herpes virus infected primary effusion lymphoma

Visit [bosterbio.com/anti-s6k1-rps6kb1-picoband-trade-antibody-a01475-2-boster.html](http://bosterbio.com/anti-s6k1-rps6kb1-picoband-trade-antibody-a01475-2-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-S6K1/RPS6KB1 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.