

Anti-CDK8 Antibody Picoband® Cy3 Conjugated

Catalog Number: A01493-1-Cy3

About CDK8

Cell division protein kinase 8 is an enzyme that in humans is encoded by the CDK8 gene. It is mapped to 13q12.13. This gene encodes a member of the cyclin-dependent protein kinase (CDK) family. CDK family members are known to be important regulators of cell cycle progression. This kinase and its regulatory subunit, cyclin C, are components of the Mediator transcriptional regulatory complex, involved in both transcriptional activation and repression by phosphorylation of the carboxy-terminal domain of the largest subunit of RNA polymerase II. This kinase regulates transcription by targeting the cyclin-dependent kinase 7 subunits of the general transcription initiation factor IIH, thus providing a link between the Mediator complex and the basal transcription machinery. Multiple pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants.

Overview

| | |
|----------------------|--|
| Product Name | Anti-CDK8 Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P49336 |

Technical Details

| | |
|------------------|---|
| Immunogen | E. coli-derived human Cdk8 recombinant protein (K355-Y464). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CDK8 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.