

Anti-Cyclin-dependent kinase 8 Cdk8 Antibody

Catalog Number: A01493-2


Overview

Product Name	Anti-Cyclin-dependent kinase 8 Cdk8 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cyclin-dependent kinase 8 Cdk8 Antibody catalog # A01493-2. Tested in WB,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49336

Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000

Anti-Cyclin-dependent kinase 8 Cdk8 Antibody (A01493-2) Images



Western blot analysis of Cdk8 on NIH/3T3 cells lysates using anti-Cdk8 antibody at 1/1,000 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cyclin-dependent kinase 8 Cdk8 Antibody

For Research Use Only. Not for use in diagnostic procedures.