

## Anti-ABCG5 Antibody

Catalog Number: A01503

### Overview

Product Name	Anti-ABCG5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ABCG5 Antibody catalog # A01503. Tested in WB, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H222

### Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 250-330
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ABCG5 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:5000 IF 1:100-300

### 1 Publications Citing This Product

1. PubMed ID: 31054325, Li R, Liu Y, Shi J, Yu Y, Lu H, Yu L, Liu Y, Zhang F. Diosgenin regulates cholesterol metabolism in hypercholesterolemic rats by

inhibiting NPC1L1 and enhancing ABCG5 and ABCG8. *Biochim Biophys Acta Mol Cell Biol Lipids*. 2019 Aug;1864(8):1124-1133. doi:10.1016/j.bbalip.2019.04.010. Epub 2019 May 1. PMID:31054325.

Visit [bosterbio.com/anti-abcg5-antibody-a01503-boster.html](https://bosterbio.com/anti-abcg5-antibody-a01503-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-ABCG5 Antibody

For Research Use Only. Not for use in diagnostic procedures.