

## Anti-Prostate Specific Antigen/KLK3 Antibody Picoband® HRP Conjugated

Catalog Number: A01505-3-HRP

### About KLK3

Prostate-specific antigen (PSA), also known as gamma-seminoprotein or kallikrein-3 (KLK3), is a glycoprotein enzyme encoded in humans by the KLK3 gene. KLK3 is a member of the kallikrein-related peptidase family and is secreted by the epithelial cells of the prostate gland. This gene is mapped to 19q13.33. KLK3 is present in small quantities in the serum of men with healthy prostates, but is often elevated in the presence of prostate cancer or other prostate disorders. KLK3 is produced for the ejaculate where it liquifies the semen in the seminal coagulum and allows sperm to swim freely. It is also believed to be instrumental in dissolving the cervical mucous cap, allowing the entry of sperm. It is not a unique indicator of prostate cancer, but may also detect prostatitis or benign prostatic hyperplasia.

### Overview

Product Name	Anti-Prostate Specific Antigen/KLK3 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P07288

### Technical Details

Immunogen	E.coli-derived human Prostate Specific Antigen/KLK3 recombinant protein (Position: A64-D255).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

## 1 Publications Citing This Product

---

1. PubMed ID: , A disposable electrochemical immunosensor based on carbon screen-printed electrodes for the detection of prostate specific antigen

Visit [bosterbio.com/anti-prostate-specific-antigen-klk3-picoband-trade-antibody-a01505-3-boster.html](https://bosterbio.com/anti-prostate-specific-antigen-klk3-picoband-trade-antibody-a01505-3-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Prostate Specific Antigen/KLK3 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.