

Anti-Cathepsin H (Y205) CTSH Antibody

Catalog Number: A01510Y205


Overview

Product Name	Anti-Cathepsin H (Y205) CTSH Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cathepsin H (Y205) CTSH Antibody catalog # A01510Y205. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P09668

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 171-210 of Human Cathepsin H.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-Cathepsin H (Y205) CTSH Antibody (A01510Y205) Images



Immunohistochemistry (IHC) analyzes of Cathepsin H (Y205) pAb in paraffin-embedded human colorectal carcinoma tissue at 1:50.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cathepsin H (Y205) CTSH Antibody

For Research Use Only. Not for use in diagnostic procedures.