

Anti-CD31/Pecam1 Antibody Picoband® Biotin Conjugated

Catalog Number: A01513-1-Biotin

About Pecam1

CD31 also known as Platelet endothelial cell adhesion molecule (PECAM-1), is a protein that in human is encoded by the PECAM1 gene. Encoded protein is a member of the immunoglobulin superfamily and this gene is mapped to 17q23.3. CD31 is found on the surface of platelets, monocytes, neutrophils, and some types of T-cells, and makes up a large portion of endothelial cell intercellular junctions. It is demonstrated that CD31 expression on human PBSCs may positively affect both neutrophil and platelet engraftment. Meanwhile, CD31 is involved in leukocyte migration and angiogenesis, which are key components of venous thrombus resolution.

Overview

Product Name	Anti-CD31/Pecam1 Antibody Picoband® Biotin Conjugated
Reactive Species	Mouse
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q08481

Technical Details

Immunogen	E. coli-derived mouse CD31 recombinant protein (Position: R58-Q273).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

20 Publications Citing This Product

1. PubMed ID: 10.4161/cbt.8.5.7687, Combined endostatin and TRAIL gene transfer suppresses human hepatocellular carcinoma growth and angiogenesis in nude mice
2. PubMed ID: 10.1155/2020/9485398, Quercetin Downregulates Cyclooxygenase-2 Expression and HIF-1alpha/VEGF Signaling-Related Angiogenesis in a Mouse Model of Abdominal Aortic Aneurysm
3. PubMed ID: -, Shen Song,Zhao Liu,Mohamed Aamer Abubaker,Ling Ding, Ji Zhang,Shengrong Yang,Zengjie Fan,Antibacterial polyvinyl alcohol/bacterial cellulose/ nano-silver hydrogels that effectively promote wound healing,Materials Science and Engineering: C,Volume 126,2021,112171,ISSN 0928-4931,https://doi.org/10.1016/j.msec.2021.112171.

Visit [bosterbio.com/anti-cd31-picoband-trade-antibody-a01513-1-boster.html](https://www.bosterbio.com/anti-cd31-picoband-trade-antibody-a01513-1-boster.html) to see all 20 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD31/Pecam1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.