

## Anti-CD31 PECAM1 Antibody

Catalog Number: A01513-4


### Overview

Product Name	Anti-CD31 PECAM1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-CD31 PECAM1 Antibody catalog # A01513-4. Tested in WB,ICC,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse.
Application	Flow Cytometry, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P16284


### Technical Details

Immunogen	This antibody is produced by immunizing rabbits with a synthetic peptide (KLH-coupled) corresponding to C-terminal residues of human CD31.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	90 mg/mL by Refractometry
Purification	Peptide affinity purified
Suggested Dilutions	WB: 1:1,000 ICC: 1:200 IHC: 1:200 FC: 1:100-1:200

## Anti-CD31 PECAM1 Antibody (A01513-4) Images



Western blot analysis of CD31 on different cell lysates using anti-CD31 antibody at 1/1000 dilution.



ICC staining CD31 in HeLa cells (green). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD31 PECAM1 Antibody

For Research Use Only. Not for use in diagnostic procedures.