

## Anti-PI3K-gamma/PIK3CG Antibody Picoband® HRP Conjugated

Catalog Number: A01517-2-HRP

### About PIK3CG

Phosphatidylinositol-4,5-bisphosphate 3-kinase catalytic subunit gamma isoform is an enzyme that in humans is encoded by the PIK3CG gene. Phosphoinositide 3-kinases (PI3Ks) phosphorylate inositol lipids and are involved in the immune response. The protein encoded by this gene is a class I catalytic subunit of PI3K. Like other class I catalytic subunits (p110-alpha p110-beta, and p110-delta), the encoded protein binds a p85 regulatory subunit to form PI3K. This gene is located in a commonly deleted segment of chromosome 7 previously identified in myeloid leukemias. Several transcript variants encoding the same protein have been found for this gene.

### Overview

Product Name	Anti-PI3K-gamma/PIK3CG Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P48736

### Technical Details

Immunogen	E.coli-derived human PI3K-gamma/PIK3CG recombinant protein (Position: K52-D308).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### 10 Publications Citing This Product

1. PubMed ID: 10.3892/mmr.2021.11938, Hydrogen sulfide ameliorates doxorubicin-induced myocardial fibrosis in rats via the PI3K/AKT/mTOR pathway
2. PubMed ID: 10.1016/j.jare.2021.09.001, A comparative study of Mesenchymal Stem Cells transplantation approach to antagonize age-associated ovarian hypofunction with consideration of safety and efficiency
3. PubMed ID: 10.1016/j.envpol.2021.116556, Mechanisms underlying reproductive toxicity induced by nickel nanoparticles identified by comprehensive gene expression analysis in GC-1 spg cells

Visit [bosterbio.com/anti-pi3k-gamma-pik3cg-picoband-trade-antibody-a01517-2-boster.html](https://bosterbio.com/anti-pi3k-gamma-pik3cg-picoband-trade-antibody-a01517-2-boster.html) to see all 10 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PI3K-gamma/PIK3CG Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.