

Anti-Toll-like receptor 8 TLR8 Antibody

Catalog Number: A01541

About TLR8

Toll-like receptors (TLRs) are signaling molecules that recognize different microbial products during infection and serve as an important link between the innate and adaptive immune responses. These proteins act through adaptor molecules such as MyD88 and TIRAP to activate various kinases and transcription factors. Like TLR7, TLR8 is localized to endosomal or lysosomal compartments and stimulates the innate immune response after activation by guanosine- and uridine-rich single-stranded RNA. Human but not murine TLR8 confers responsiveness to the antiviral compound R-848.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-Toll-like receptor 8 TLR8 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Toll-like receptor 8 TLR8 Antibody (Catalog # A01541). Tested in ELISA, IHC-P applications. This antibody reacts with Human, Mouse.
Application	ELISA, IHC-P
Clonality	Polyclonal
Formulation	TLR8 Antibody is supplied in PBS containing 0.02% sodium azide.
Storage Instructions	TLR8 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. Avoid repeated freeze-thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Host	Rabbit
Uniprot ID	Q9NR97

Technical Details


Immunogen	TLR8 antibody was raised against a 13 amino acid synthetic peptide near the amino terminus of human TLR8. The immunogen is located within the first 50 amino acids of TLR8.
Predicted Reactive Species	Bovine
Cross Reactivity	At least three isoforms of TLR8 are known to exist. TLR8 antibody is predicted to not cross-react

	with other TLR protein family members
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL
Purification	TLR8 Antibody is affinity chromatography purified via peptide column.
Suggested Dilutions	TLR8 antibody can be used for detection of TLR8 by immunohistochemistry at 5 ug/mL. Antibody can also be used for immunohistochemistry starting at 5 ug/mL. Antibody validated: Immunohistochemistry in human and mouse samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.

Anti-Toll-like receptor 8 TLR8 Antibody (A01541) Images



Immunohistochemistry of TLR8 in human thymus tissue with TLR8 antibody at 5 ug/mL.



Immunohistochemistry of TLR8 in mouse thymus tissue with TLR8 Antibody at 5 ug/mL.

[Anti-Toll-like receptor 8 TLR8 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.