

## Anti-Ribosomal Protein S3 RPS3 Antibody

Catalog Number: A01542-1

### Overview

Product Name	Anti-Ribosomal Protein S3 RPS3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Ribosomal Protein S3 RPS3 Antibody catalog # A01542-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal 3A1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P23396

### Technical Details

Immunogen	Synthesized peptide derived from the C-terminal region of human Ribosomal Protein S3. at AA range: 140-220
Predicted Reactive Species	Bovine, Canine, Equine, Guinea Pig, Rabbit, Sheep, Zebrafish
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	This antibody's concentration is >0.5mg/ml.
Purification	Immunogen affinity purified
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB, 1:500-1:2000</p>

## Anti-Ribosomal Protein S3 RPS3 Antibody (A01542-1) Images



Western Blot (WB) analysis of specific cells using Ribosomal Protein S3 Polyclonal antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ribosomal Protein S3 RPS3 Antibody