

Anti-Fos B (S23) Antibody

Catalog Number: A01569S23


Overview

Product Name	Anti-Fos B (S23) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Fos B (S23) Antibody catalog # A01569S23. Tested in WB,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P53539


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids N-terminus of Human Fos B.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IF: 1:50-1:200

Anti-Fos B (S23) Antibody (A01569S23) Images



Western blot (WB) analysis of Fos B (S23) polyclonal antibody at 1:500 dilution Lane1:HEK293T whole cell lysate Lane2:sp2/0 whole cell lysate Lane3:PC12 whole cell lysate



Immunohistochemistry (IHC) analyzes of FosB (S23) pAb in paraffin-embedded human prostate carcinoma tissue.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fos B (S23) Antibody

For Research Use Only. Not for use in diagnostic procedures.