

## Anti-Caveolin-2 (phospho-Y19) CAV2 Antibody

Catalog Number: A01574Y19

### Overview


|                      |  |
|----------------------|--|
| Product Name         | Anti-Caveolin-2 (phospho-Y19) CAV2 Antibody  |
| Reactive Species     | Human  |
| Description          | Boster Bio Anti-Caveolin-2 (phospho-Y19) CAV2 Antibody catalog # A01574Y19. Tested in IHC applications. This antibody reacts with Human.   |
| Application          | IHC  |
| Clonality            | Polyclonal   |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2   |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | P51636   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Synthetic phosphopeptide derived from human Caveolin-2 around the phosphorylation site of Tyrosine19.   |
| Cross Reactivity    | No cross reactivity with other proteins.  |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | IHC: 1:50-1:200   |

## Anti-Caveolin-2 (phospho-Y19) CAV2 Antibody (A01574Y19) Images

---



Immunohistochemistry (IHC) analyzes of Caveolin-2 (phospho-Y19) polyclonal antibody in paraffin-embedded human prostate tissue at 1:50.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Caveolin-2 (phospho-Y19) CAV2 Antibody

For Research Use Only. Not for use in diagnostic procedures.