

## Anti-Rad52 (phospho-Y104) Antibody

Catalog Number: A01580Y104


### Overview

Product Name	Anti-Rad52 (phospho-Y104) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Rad52 (phospho-Y104) Antibody catalog # A01580Y104. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P43351

### Technical Details

Immunogen	Synthetic phosphopeptide derived from human Rad52 around the phosphorylation site of Tyrosine 104.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-Rad52 (phospho-Y104) Antibody (A01580Y104) Images



Western blot (WB) analysis of p-Rad52 (Y104) polyclonal antibody at 1:500 dilution  
Lane1:HepG2 cell lysate  
Lane2:Mouse spleen tissue lysate Lane3:Rat spleen tissue lysate

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rad52 (phospho-Y104) Antibody

For Research Use Only. Not for use in diagnostic procedures.