

## Anti-Claudin 1/CLDN1 Antibody Picoband® HRP Conjugated

Catalog Number: A01585-3-HRP

### About CLDN1

Claudin-1 is a protein that in humans is encoded by the CLDN1 gene. Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell sheets, forming continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. These junctions are comprised of sets of continuous networking strands in the outwardly facing cytoplasmic leaflet, with complementary grooves in the inwardly facing extracytoplasmic leaflet. The protein encoded by this gene, a member of the claudin family, is an integral membrane protein and a component of tight junction strands. Loss of function mutations result in neonatal ichthyosis-sclerosing cholangitis syndrome.

### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-Claudin 1/CLDN1 Antibody Picoband® HRP Conjugated                              |
| Reactive Species     | Human   |
| Clonality            | Polyclonal  |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.    |
| Host                 | Rabbit  |
| Uniprot ID           | O95832  |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | E.coli-derived human Claudin 1/CLDN1 recombinant protein (Position: F96-V211).                     |
| Cross Reactivity    | No cross-reactivity with other proteins.   |
| Isotype             | Rabbit IgG   |
| Form                | Liquid   |
| Concentration       | 0.5 mg/mL  |
| Purification        | Immunogen affinity purified.   |
| Conjugate           | HRP  |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Claudin 1/CLDN1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.