

Anti-Parathyroid Hormone Receptor 1/PTH1R Antibody Picoband® Biotin Conjugated

Catalog Number: A01588-Biotin

About PTH1R

Parathyroid hormone/parathyroid hormone-related peptide receptor, also known as PTH1R, PTHR, is a protein that in humans is encoded by the PTH1R gene. PTHR is a member of a family of G protein-coupled receptors. PTHR gene was mapped to the vicinity of the 3p21.3-p21.2 boundary by PCR analysis of human/rodent somatic cell hybrid panels using oligonucleotide primers designed to amplify a portion of the gene from genomic DNA. This is a receptor for parathyroid hormone and for parathyroid hormone-related peptide. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase and also a phosphatidylinositol-calcium second messenger system. PTH1R functions as a receptor for parathyroid hormone and for parathyroid hormone-related protein.

Overview

Product Name	Anti-Parathyroid Hormone Receptor 1/PTH1R Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Monkey
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q03431

Technical Details

Immunogen	E.coli-derived human Parathyroid Hormone Receptor 1/PTH1R recombinant protein (Position: R102-L446).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental

requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Parathyroid Hormone Receptor 1/PTH1R Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.