

## Anti-NEDD4L (phospho-S448) Antibody

Catalog Number: A01595S448


### Overview

Product Name	Anti-NEDD4L (phospho-S448) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NEDD4L (phospho-S448) Antibody catalog # A01595S448. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96PU5

### Technical Details

Immunogen	A synthetic phosphorylated peptide around S448 of human NEDD4L
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

## Anti-NEDD4L (phospho-S448) Antibody (A01595S448) Images



Western blot analysis of extracts of various cell lines, using Phospho-NEDD4L-S448 antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NEDD4L (phospho-S448) Antibody

For Research Use Only. Not for use in diagnostic procedures.