

Anti-eNOS (Phospho-Ser1179) NOS3 Antibody

Catalog Number: A01604S1179

Overview

Product Name	Anti-eNOS (Phospho-Ser1179) NOS3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-eNOS (Phospho-Ser1179) NOS3 Antibody catalog # A01604S1179. Tested in IHC, IF, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P29474

Technical Details

Immunogen	Synthesized peptide derived from human eNOS (Phospho-Ser1179)
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	IHC 1:50-200 WB 1:500-2000 IF 1:50-200

1 Publications Citing This Product

1. PubMed ID: 33933809, Zhang J,Cui J,Li X,Hao X,Guo L,Wang H,Liu H.Increased secretion of VEGF-C from SiO₂-induced pulmonary macrophages promotes lymphangiogenesis through the Src/eNOS pathway in silicosis.Ecotoxicol Environ Saf.2021 Apr 29;218:112257.doi:10.1016/j.ecoenv.2021.112257.Epub ahead of print.PMID:33933809.

Visit bosterbio.com/anti-enos-phospho-ser1179-nos3-antibody-a01604s1179-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-eNOS (Phospho-Ser1179) NOS3 Antibody

For Research Use Only. Not for use in diagnostic procedures.