

## Anti-NEK2 Antibody

Catalog Number: A01606-1


### Overview

Product Name	Anti-NEK2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-NEK2 Antibody catalog # A01606-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P51955


### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence within amino acids 345-445 of human NEK2 (NP_002488.1).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:2000-1:5000 IHC: 1:50-1:200

## Anti-NEK2 Antibody (A01606-1) Images



Western blot analysis of lysates from HeLa cells, using NEK2 pAb at 1:1000 dilution.



Immunohistochemistry analysis of NEK2 in paraffin-embedded rat ovary using NEK2 pAb at dilution of 1:100 (40x lens).

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-NEK2 Antibody

For Research Use Only. Not for use in diagnostic procedures.