

## Anti-SLK Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A01614-2-Fluoro550

### About SLK

STE20-like serine/threonine-protein kinase is an enzyme that in humans is encoded by the SLK gene. This gene is mapped to chromosome 10q24.33-q25.1. SLK is a kinase with multiple functional domains that is predicted to regulate cytoskeletal organization and responses to apoptotic stimuli.

### Overview

Product Name	Anti-SLK Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9H2G2

### Technical Details

Immunogen	E.coli-derived human SLK recombinant protein (Position: E328-S1232). Human SLK shares 78.2% and 72.1% amino acid (aa) sequence identity with mouse and rat SLK, respectively.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLK Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.