

Anti-Mucin-16 MUC16 Antibody

Catalog Number: A01641


Overview

Product Name	Anti-Mucin-16 MUC16 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Mucin-16 MUC16 Antibody catalog # A01641. Tested in WB,ICC/IF,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8WXI7


Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 FC: 1:50-1:100 ICC: 1:100-1:500

Anti-Mucin-16 MUC16 Antibody (A01641) Images



ICC staining MUC16 in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining MUC16 in Hela cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

1 Publications Citing This Product

1. PubMed ID: 32765681, Wang H, Yan C, Ye H. Overexpression of MUC16 predicts favourable prognosis in MUC16-mutant cervical cancer related to immune response. *Exp Ther Med*. 2020 Aug;20(2):1725-1733. doi:10.3892/etm.2020.8836. Epub 2020 Jun 4. PMID:32765681; PMCID:PMC7388571.

Visit bosterbio.com/anti-mucin-16-muc16-antibody-a01641-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mucin-16 MUC16 Antibody

For Research Use Only. Not for use in diagnostic procedures.