

## Anti-Tryptophan 5-hydroxylase 2 TPH2 Antibody

Catalog Number: A01643

### Overview


Product Name	Anti-Tryptophan 5-hydroxylase 2 TPH2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Tryptophan 5-hydroxylase 2 TPH2 Antibody catalog # A01643. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8IWU9

### Technical Details

Immunogen	Recombinant full length Human TPH2.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200 IF: 1:50-1:200

## Anti-Tryptophan 5-hydroxylase 2 TPH2 Antibody (A01643) Images

---



Immunofluorescence analysis of A549 cell using TPH2 antibody. Blue: DAPI for nuclear staining.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Tryptophan 5-hydroxylase 2 TPH2 Antibody

For Research Use Only. Not for use in diagnostic procedures.